Nutri-Bact

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BHI Infusion Broth with 6,5% Sodium Chloride is a liquid medium used to differentiate enterococci from other streptococci by addition of 68 grams per liter of sodium chloride. Sodium chloride reacts as selective agents because mostly of Enterococci sp. may grow in this selective broth. Bromocresol purple is the indicator.

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This medium is for *in vitro* diagnostic use only.

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Store prepared media between 2-8°C, protected from direct light. Store dehydrated media in a dry place, in its tightly-sealed original container between 2-25°C.

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Media should not be used if the expiry date has passed. Prepared media should not be used if there are signs of contamination or deterioration (shrinking, cracking, evaporation or discoloration). Do not use dehydrated media if it is caked.

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Suspend **91,0 g** of dehydrated medium in a final volume of 1000 ml distilled water. Heat gently with frequent agitation to dissolve completely. Distribute into culture tubes and autoclave at 121°C for 15 minutes.

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Prior to inoculate, the prepared media should be brought to room temperature.

1. Inoculate with 3-5 similar colonies from an overnight growth.
2. Incubate at 35°C for 48-72 hours.
3. Examine daily for growth.

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Results after 24 hours at  $35 \pm 2^\circ\text{C}$  aerobically.

Souches / Strains	ATCC	Croissance / Growth	Couleur / Color
<i>Enterococcus faecalis</i>	29212	+	Rouge / Red
<i>Streptococcus pyogenes</i>	19615	-	Bleu / blue
<i>Streptococcus bovis</i>	14028	-	Bleu / blue

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This medium allows partial identification. Additional biochemical tests, as well as morphological characteristics and serological typing, may be required. Negative result must be confirmed only after 3 days of incubation.

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1. Facklam, R.R., 1973 Comparaison of several laboratory media for presumptive identification of enterococci and group D streptococci. Appl. Microbiol. 26:138.
2. Lennette, E.H., Ballows, A., Hausler, W.J.Jr., and Shadomy, H.J. Manual of Clinical Microbiology. 4th ed. 1985 Washington D.C.: American society for Microbiology.
3. Mac Faddin, Jean F. 1985 Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Vol. I, Baltimore, MD.: Williams & Wilkins.

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