Nutri-Bact

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Bile Esculin Agar is a selective culture medium used for the isolation of group D streptococci. Esculin hydrolysed by group D gives a dark brown colour. Oxgall inhibits Gram-positive bacteria other than enterococci.

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This medium is for *in vitro* diagnostic use only.

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Store prepared media between 2-8°C, protected from direct light. Store dehydrated media in a dry place, in its tightly-sealed original container between 2-25°C.

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Media should not be used if the expiry date has passed. Prepared media should not be used if there are signs of contamination or deterioration (shrinking, cracking, evaporation or discoloration). Do not use dehydrated media if it is caked.

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Suspend **64.0 g** of **Bile Esculin Agar** powder into a final volume of 1000 ml of filtered deionized water. Heat with frequent agitation in order to dissolve completely, then let boil for one minute. Autoclave at 121°C, 15 p.s.i. during 15 minutes. Once sterilized, cool down to 45-50°C. Add aseptically 50 ml of sterile horse serum (#4173) if wanted. Mix gently and dispense aseptically.

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Prior to inoculate, the prepared media should be brought to room temperature.

1. Inoculate with 2-3 similar colonies.
2. Incubate at 35°C for 24 hours.

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1. Lennette, E.H., Ballows, A., Hausler, W.J.Jr., and Shadomy, H.J. Manual of Clinical Microbiology. 4th Ed. 1985 Washington D.C.: American society for Microbiology.
2. Mac Faddin, Jean F. 1985 Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Vol. I, Baltimore, MD.: Williams & Wilkins.
3. N.C.C.L.S. 1990 Quality Assurance for Commercially Prepared Microbiological Culture Media Approved Standard. Vol. 10, No. 14, NCCLS Document M22-A.

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