Nutri-Bact

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Hektoen Agar is a selective and differential solid medium for the cultivation and isolation of Gram-negative enteric microorganisms from clinical and nonclinical specimens. Addition of sucrose, salicin and lactose differentiate enteric pathogens from slow-lactose-fermenting organisms. Bile salts inhibit gram-positive. Complex sodium thiosulfate-ferric ammonium citrate produces black-centered colonies when ferric ammonium citrate precipitates (H<sub>2</sub>S production).

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This medium is for *in vitro* diagnostic use only.

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Store prepared media between 2-8°C, protected from direct light. Store dehydrated media in a dry place, in its tightly-sealed original container between 2-25°C.

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Media should not be used if the expiry date has passed. Prepared media should not be used if there are signs of contamination or deterioration (shrinking, cracking, evaporation or discoloration). Do not use dehydrated media if it is caked.

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Suspend **75 g** of dehydrated medium in a final volume of 1000 ml distilled water. Heat gently with frequent agitation to dissolve completely and boil for one minute. DO NOT AUTOCLAVE. Cool to 45-50° C. Mix gently and dispense aseptically.

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Prior to inoculate, the prepared media should be brought to room temperature.

1. Directly, inoculate plate as soon as possible with specimen and streak for isolation. If inoculation will be delayed, an Enteric Pathogen Transport Medium as EPT-100 should be used.
2. Incubate at  $35 \pm 2^\circ\text{C}$  for 24-48 hours.

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Results after 24 hours at  $35 \pm 2^\circ\text{C}$  aerobically.

Souches / Strains	ATCC	Croissance / Growth	Couleur / Color
<i>Salmonella typhimurium</i>	14028	+	Bleu + H <sub>2</sub> S / Blue + H <sub>2</sub> S
<i>Shigella flexneri</i>	12022	+	Bleu + H <sub>2</sub> S / Blue + H <sub>2</sub> S
<i>Escherichia coli</i>	25922	- ou partielle / or partial	
<i>Streptococcus faecalis</i>	29212	- ou partielle / or partial	

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This medium allows partial identification. Additional biochemical tests, as well as morphological characteristics and serological typing, may be required.

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1. Lennette, E.H., Ballows, A., Hausler, W.J.Jr. and Shadomy, H.J. Manual of Clinical Microbiology. 4th ed. 1985 Washington D.C.: American society for Microbiology.
2. N.C.C.L.S. 1990 Quality Assurance for Commercially Prepared Microbiological Culture Media. Approved Standard. Vol.10. No.14. NCCLS Document M22-A.

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3. Mac Faddin, Jean F., 1985 Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria. Vol.1. Baltimore, MD.: Williams & Wilkins.

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