Nutri-Bact

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Sabouraud Dextrose Agar is a solid medium used for isolation, cultivation and maintenance of pathogenic and non-pathogenic yeasts and molds. The high dextrose concentration and acidic pH permits selectivity of fungi.

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This medium is for *in vitro* diagnostic use only.

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Store prepared media between 2-8°C, protected from direct light. Store dehydrated media in a dry place, in its tightly-sealed original container between 2-25°C.

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Media should not be used if the expiry date has passed. Prepared media should not be used if there are signs of contamination or deterioration (shrinking, cracking, evaporation or discoloration). Do not use dehydrated media if it is caked.

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Suspend **67.0 g** into a final volume of 1000 ml of filtered deionized water. Heat with frequent agitation in order to dissolve completely, then let boil for one minute. Autoclave at 121°C during 15 minutes. Once sterilized, cool down to 45-50°C. Mix gently and dispense aseptically.

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Prior to inoculate, the prepared media should be brought to room temperature.

1. Directly, inoculate plates and streak for isolated colonies. Inoculate each specimen in duplicate.
2. Seal plates with masking tape, incubate a set at 35 ± 2°C for 7 days and the other set, at room temperature for up to 21 days. Tubes should be incubated caps loosened.
3. Examine frequently every day or every 2-3 days for growth.

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Results after 7 days at 25°C aerobically.

Souches / Strains	ATCC	Croissance / Growth
<i>Candida albicans</i>	60193	+
<i>Candida albicans</i>	10231	+
<i>Saccharomyces cerevisiae</i>	9763	+
<i>Aspergillus brasiliensis</i>	16404	+
<i>Trichophyton mentagrophytes</i>	9533	+

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This medium allows partial identification. Additional biochemical tests, as well as morphological characteristics and serological typing, may be required.

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1. Lennette, E.H., Ballows, A., Hausler, W.J. Jr, and Shadomy, H.J. 1985. Manual of Clinical Microbiology, 4th Ed., Washington, DC. American Society for Microbiology.
2. Mac Faddin, Jean, F. 1985. Identification-Cultivation-Maintenance of Medical Bacteria Vol.1, Baltimore, MD, Williams & Wilkins.
3. N.C.C.L.S. 1996. Quality Assurance for Commercially Prepared Microbiological Culture Media Approved Standard. N.C.C.L.S. Document M22-A2 Vol.16, No. 16.

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