Nutri-Bact

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OF Medium is a semisolid medium, supplemented with carbohydrate, for determination of the oxidative-fermentative metabolism of this carbohydrate by gram negative bacilli. Motility and gas-production may be determined due to the semisolid nature of the medium. Tubes should be inoculated in duplicate: cover one with 1-2 ml of sterile mineral oil for the fermentative metabolism.

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This medium is for *in vitro* diagnostic use only.

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Store prepared media between 2-8°C, protected from direct light. Store dehydrated media in a dry place, in its tightly-sealed original container between 2-25°C.

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Media should not be used if the expiry date has passed. Prepared media should not be used if there are signs of contamination or deterioration (shrinking, cracking, evaporation or discoloration). Do not use dehydrated media if it is caked.

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Suspend **9,8 g** of dehydrated media in 900 ml of purified filtered water. Heat with frequent agitation and boil for one minute. To prepare OF 1% of the desired carbohydrate, add 100 ml of 10% carbohydrate solution. To prepare OF Control, add 100 ml of purified filtered water. Sterilize at 121° C for 15 minutes. Cool to 45-50° C. Mix gently and dispense into sterile culture tubes. Cool in upright position.

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Prior to inoculate, the prepared media should be brought to room temperature.

1. From a pure overnight culture, prepare a bacterial heavy suspension.
2. Inoculate 2 tubes with a straight needle not far than 1/2 cm from the bottom of the tubes.
3. Add 1-2 ml of sterile mineral oil to one of the 2 tubes. Loosen cap of the other tube
4. Incubate at 35° C for 24-48 hours. Discard only after 4 days of incubation

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This medium allows partial identification. Additional biochemical tests, as well as morphological characteristics and serological typing, may be required.

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1. Hugh, R., E. Leifson 1953 The taxonomic significance of fermentative versus oxidative metabolism of carbohydrates by various gram-negative bacteria. J. Bacteriol. 66:24-26
2. Lennette, E.H., Ballows, A., Hausler, W.J.Jr., and Shadomy, H.J. Manual of Clinical Microbiology. 4th ed. 1985 Washington D.C.: American society for Microbiology.
3. Mac Faddin, Jean F., Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria Vol.1 1985 Baltimore, MD. Williams & Wilkins.
4. Washington, J.A. 1981 Laboratory Procedures in Clinical Microbiology Springer-Verlag, New York

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