

STOOL SPECIMEN COLLECTION AND TRANSPORT SYSTEMS

EPT-120

ENTERIC PATHOGEN
TRANSPORT MEDIUM

SAF-120

FIXATIVE FOR PARASITOLOGY
(SODIUM ACETATE + ACETIC ACID + FORMALINE)



POLYPROPYLENE VIAL OF 30 ML CONTAINING 15 ML OF EPT MEDIUM OR SAF FIXATIVE OR OTHER, WITH A COLORED SCREW CAP AND LABEL CLEARLY MARKED "Fill up to line", CONTAINING A SCOOP BUILT INTO THE CAP DESIGNED TO FACILITATE THE COLLECTION AND HANDLING OF SAMPLES FROM PATIENTS TO INCREASE THE ACCURACY OF DIAGNOSIS.

DESIGNED FOR THE SAFETY OF LABORATORY AND NURSING PERSONNEL, PARTICULAR CARE WAS GIVEN TO ELIMINATE LEAKAGES DURING TRANSPORT BY APPLYING STRICT QUALITY CONTROL AND INSPECTION STEPS DURING THE MANUFACTURING PROCESS.

NUTRI-BACT OFFERS THE TWO MOST COMMONLY USED COLLECTION AND TRANSPORT SYSTEMS, EACH AVAILABLE IN BOX OF **25 OR 120 UNITS.**